

```

hs_sch 'lt5'; # we choose the sample changer lt5
# everything after #-symbol is a comment
# a3 is at -15

hsarot 'on'; # sample rotation is ON,
ttol 3;      # default temperature tolerance

# technical code for the temperature regulation
# with new double-1-active heater
tdbl 1;
sics '
table fix_warmup_weight 1
warmup weight -0.6
warmup trig 10
warmup limit 70
warmup timef 0.2
warmup abruptstop 1
';
# END of technical code for the temperature regulation

hs_lambda '1p15'; hs_resol 'MR';
# we choose 1.15A and Medium Resolution (MR)
# hs_lambda '1p89'; hs_resol 'HI';
# possible choice with 1.89A and
# High Intensity (HI)

psinq 1450; # really optional to readjust the presets
# for SING current 1450uA. The default table values are
# for 1500uA. This adjustment tunes one sweep to be
# closer to 1/2h. Might be important if SING
# current is very different from 1500uA.

# psweep 1; # time in hours for one sweep. Default value is
# 1/2 hour.

for ($t=230; $t<300.1; $t += 15) {
# this is the loop over the temperatures

hcount 6, 'La0.7Sr0.3MnO3_2g,V6x20',1,$t,'600 3';
# we count 6 sweeps (default sweep time is 1/2h)
# with the sample in position no.1 at
# the temperatures 230,245, ... After the temperature
# is in tolerance 3K with the set-point we wait 600s and
# then count.
}

hcount 300, 'La0.7Sr0.3MnO3_2g,V6x20', 1, 300, '600 5';

# this is the end of the script.
__END__
The __END__ is optional, but can be useful...
Everything after the above __END__ is ignored. Might be used for some
chunks of code you would like to keep, avoiding commenting each line
with '#'-symbol

# hs_resol 'MR';
# hs_lambda '1p89';
# hcount 10, 'La0.7Sr0.3MnO3_2g,V6x20',1,1;

```